

Student: JAYSON YATES

Aseptic transfer technique for agar slant: correct order of steps

1. Pick up and heat loop
2. Pick up culture, remove cap, heat mouth of tube
3. Insert loop into culture
4. Heat mouth of culture tube, replace cap, place in rack
5. Pick up sterile slant, remove cap, heat mouth of tube
6. Streak culture onto sterile slant with loop
7. Heat mouth of sterile slant, replace cap, place in rack
8. Heat loop and return to rack

Collected Lab Data

Tube	Negative control	<i>E. coli</i> inoculation
Growth result	No growth	Growth
Image		